Cooltouch II Specifications

1320nm Wavelength

- Designed specifically for its absorption in water
- Allows for deep penetration and uniform treatment
- Easily treats all skin types without pigmentary changes

Thermal Feedback

- Measures immediate response to ensure therapeutic treatment levels are achieved
- Serves as a gauge to vary the laser energy delivered to different areas
- Allows you to customize treatments for each patient for desired results

Precision Cooling

- 2 settings: Variable pre and post-cooling modes maximize patient comfort
- Soothes and protects the skin's surface during treatment

Voltage: 110V-AC or 220V-AC